



Razorcat Development GmbH
Witzlebenplatz 4
14057 Berlin
Germany

Technical Support

phone	+49 - (30) - 536 357 0
fax	+49 - (30) - 536 357 60
email	support@razorcat.com
www	http://www.razorcat.com

Release Notes

Release v4.3.1 (15 items)

CTE Common

New Feature

CTE-84

Added option to export CTE documents as PNG, SVG and JPG images.

CTE Common

New Feature

CTE-93

Added search and replace dialog to find strings in CTE document components.

CTE Common

Enhancement

CTE-375

In exported diagrams, connection lines are now displayed behind the nodes. When editing, they remain in the foreground for more intuitive editing.

CTE Common

New Feature

CTE-455

Allow copy and paste between the standalone CTE and the CTE integrated in TESSY 5.1.9 and later.

CTE Common

New Feature

CTE-570

Added option to export CTE documents to a stable XML format for processing in other tools.

Combination Mark Editing

Enhancement

CTE-152

Improved selection in the Combination Table, e.g. "Select All" (CTRL+A) now selects all marks in the same test table column as the currently selected mark.

Export

New Feature

CTE-529

Added option to export CTE documents in an EXCEL tree structure.

Release v4.3.1 (15 items)

Layout

Enhancement

CTE-332

The Layout button now has an additional menu item that triggers a layout of the entire document, i.e. all open and closed refinements of the document are laid out.

Outline

New Feature

CTE-287

The refinement hierarchy of a CTE document is now also displayed in the outline view.

Table Editing

New Feature

CTE-149

Tree nodes are highlighted when the mouse is over a mark. Tree nodes associated with the selection are also highlighted.

Table Editing

New Feature

CTE-171

It is now possible to show/hide/create/delete/move automatic and manual page breaks in the test table to predefine the content of the table pages. The image export creates a separate image for each table page.

Table Editing

Enhancement

CTE-414

CTE now synchronizes test case selection across all refinements.

Tree Editing

Enhancement

CTE-54

The user can now define font styles (type, weight, size) and text color for tree nodes and table entries.

Tree Editing

Enhancement

CTE-276

Display of nodes that are not connected to the root node under a virtual node "Nodes not in Tree".

Tree Editing

Enhancement

CTE-582

Added automatic switching between quick edit mode and tree editing. This reduces the need to switch from mouse and back when creating and editing trees.